


A Brachial Plexus Contouring Guide

DOSEPEDIA

1

DOSEPEDIA

Introduction



2


Introduction Cont.

- The Brachial Plexus can be difficult to contour on a CT scan.
- At the time of this presentation, there currently lacks a thorough guide with visual aids to contour the Brachial Plexus intended for Physics staff.
- This guide will involve three distinct steps which when completed and used in combination, will allow for a reliable contour of the Brachial Plexus.
 1. Contouring the Anterior and Middle Scalene Muscles.
 2. Using the Timmerman Technique.
 3. Following RTOG's 'Brachial Plexus Contouring Atlas'.
- However, if you can see it on the CT, then contour it!
- Slice-by-slice contouring can be found at the end of the presentation.

3

Brachial Plexus Contour Overview

- The plexus extends from the spinal cord and is comprised of nerves: C5, C6, C7, C8, & T1.
- On the frontal view, the nerves appear to extend from in between the vertebral bodies (VB):
 - VBC4 - C5 - VBC5
 - VBC5 - C6 - VBC6
 - VBC6 - C7 - VBC7
 - VBC7 - C8 - VBT1
 - VBT1 - T1 - VBT2
- RTOG's 'Brachial Plexus Contouring Atlas': 4 nerves.




4

DOSEPEDIA

Step 1: Contouring the Scalene Muscles


- Anterior scalene muscle (ASM):
 - Inserts into the scalene tubercle of the first rib.
 - Arises from the anterior tubercles of the transverse processes of C3 to C6.
- Middle scalene muscle (MSM):
 - Inserts into the upper surface of the first rib.
 - Arises from the posterior tubercles of the transverse processes of C2 to C7
- **Use the frontal view to contour.**
 - Reference the axial view to tweak contours.



5

DOSEPEDIA


Step 1: Anterior Scalene Muscle



6

DOSEPEDIA


Step 1: Anterior Scalene Muscle Cont.



7

DOSEPEDIA


Step 1: Anterior Scalene Muscle Cont.



8

DOSEPEDIA


Step 1: Anterior Scalene Muscle Cont.



9

DOSEPEDIA


Step 1: Anterior Scalene Muscle Cont.



10

DOSEPEDIA


Step 1: Anterior Scalene Muscle Cont.



11

DOSEPEDIA


Step 1: Middle Scalene Muscle



12

DOSEPEDIA


Step 1: Middle Scalene Muscle Cont.



13

DOSEPEDIA


Step 1: Middle Scalene Muscle Cont.



14

DOSEPEDIA


Step 1: Middle Scalene Muscle Cont.



15

DOSEPEDIA


Step 1: Middle Scalene Muscle Cont.



16

DOSEPEDIA

Step 1: Middle Scalene Muscle Cont.



17

DOSEPEDIA

Step 2: Using the Timmerman Technique

Locating the Brachial Plexus
Timmerman's Trick-1

A coronal CT scan image of the upper thorax and neck. The brachial plexus is highlighted with yellow contours. Labels 'clavicle' and '1st rib' point to the respective anatomical structures. The yellow contours represent the rough outline of the neurovascular bundle as described in the text.

- Vein, artery, and nerve (VAN, anterior to posterior) will go over the 1st rib and under the clavicle
- Using coronal images, find the plane where vascular/nerve structures (tubes and wires) pass between the 1st rib and clavicle
- Roughly contour these neurovascular tissues in this coronal plane (as shown in yellow)
- You will use these rough contours in the next step


RTOG 20

18

DOSEPEDIA

Step 2: Using the Timmerman Technique Cont.


- Approximates for the path of the nerve roots.
- Use the frontal view to contour.
 - Parallel to the vertebral column.
 - Along the arc of the chest wall.
- Using a static 3D-Brush, contour on the slice where both scalene muscles are present.
- Reference the axial view to verify that the contour is reasonable.



19

DOSEPEDIA


Step 2: Using the Timmerman Technique Cont.



20

DOSEPEDIA


Step 2: Using the Timmerman Technique Cont.



21

DOSEPEDIA


Step 1 & 2: ASM, MSM, & Timmerman Tech.



22

DOSEPEDIA

Step 1 & 2: ASM, MSM, & Timmerman Tech. Cont.




23

DOSEPEDIA

Step 3: Contouring via RTOG Guidelines

“ ...

- Identify and contour C5, T1, and T2.
- Identify and contour the subclavian and axillary neurovascular bundle.
- Identify and contour anterior and middle scalene muscles from C5 to insertion onto the first rib.
- To contour the brachial plexus OAR use a 5-mm diameter paint tool.
- Start at the neural foramina from C5 to T1; this (BP) should extend from the lateral aspect of the spinal canal to the small space between the anterior and middle scalene muscles.




24

DOSEPEDIA

Step 3: Contouring via RTOG Guidelines Cont.

- For CT slices, where no neural foramen is present, contour only the space between the anterior and middle scalene muscles.
- Continue to contour the space between the anterior and middle scalene muscles; eventually the middle scalene will end in the region of the subclavian neurovascular bundle.
- Contour the brachial plexus as the posterior aspect of the neurovascular bundle inferiorly and laterally to one to two CT slices below the clavicular head.
- The first and second ribs serve as the medial limit of the OAR contour.

“ ... ”



25


DOSEPEDIA

Step 3: Contouring via RTOG Guidelines

“ ... ”

- Identify and contour C5, T1, and T2

“ ... ”



26


DOSEPEDIA

Step 3: Contouring via RTOG Guidelines Cont.

“ ...

- Start at the neural foramina, and extend from the spinal canal to the space between the ASM & MSM

” ...



27


DOSEPEDIA

Step 3: Contouring via RTOG Guidelines Cont.

“ ...

- Where no neural foramen is present, contour only the space between the ASM & MSM

” ...



28

Eclipse Shortcut Keys, & Contouring Tips & Tricks

DOSEPEDIA

29

DOSEPEDIA

Contouring Shortcut Keys

- Rotate through different plane views.
 - 'control' key & 'r'.



- 'control' key - pinky
- 'r' - index



30


DOSEPEDIA

Contouring Shortcut Keys Cont.

- Hide all contours other than the highlighted structure.
 - Hold 'h'.



- 'h' - index



31


DOSEPEDIA

Contouring Shortcut Keys Cont.

- 'control' key + 'shift' key + **1/2/6/7**.
- Compare top right screens in External Beam Planning in Plan Evaluation (2).
- View a plan as if in External Beam Planning in Plan Evaluation (6).
- Present up to 6 plans and their respective DVHs simultaneously (7).



- 'control' key - pinky
- 'shift' key - ring



32


DOSEPEDIA

Contouring Shortcut Keys Cont.

- Reload all.
 - 'alt' key → 'f' → 'r'.



- 'alt' key - thumb
- 'f' - index
- 'r' - middle



33

DOSEPEDIA

Contouring Shortcut Keys Cont.

- Scroll through CT while skipping slices.
 - Hold 'alt' key & Scroll 'mouse wheel'.



- 'alt' key - thumb
- 'scroll' - right middle



34

Contouring Shortcut Keys Cont.

- Following an image fusion, blend fusion to favor the CT or fused image without the slider.
 - Hold 'control' key & 'a'.
 - The image blend must be initially set to a percentage other than 50%, and doing so will change the current blend to the complement of the initial blend percentage.



- 'control' key - pinky
- 'a' - ring or middle

35

Contouring Shortcut Keys Cont.

- Create a new structure.
 - Display the screen to create a new structure.
 - 'alt' key → 's' → 'enter' key.
 - Fill in structure information.
 - 'down' key → 'tab' key → Type Structure Name → 'enter' key.



- Display window level.
 - 'alt' key → 'v' → 'enter' key.



- 'alt' key - thumb
- 's' - ring
- 'down' key - pinky
- 'tab' key - right ring or pinky

- 'alt' key - thumb
- 'v' - index
- 'enter' key - right thumb

36

DOSEPEDIA

Contouring Shortcut Keys Cont.

- Change window level without the slider by using only the mouse cursor.
 - Select Mouse Cursor Tool → Hold 'shift' key & Move Mouse Cursor.
 - Left/Right: Tighten/widen contrast boundaries.
 - Up/Down: Translate contrast boundaries.
 - Diagonal: Combination.



- 'shift' key - pinky or ring

- Auto window level adjustment, and useful as a contrast reset.
 - 'alt' key → 'v' → 'a'.



- 'alt' key - thumb
- 'v' - index
- 'a' - pinky

37

DOSEPEDIA

Contouring Shortcut Keys Cont.




38

Adaptive Brush Contouring

- This method involves manipulating the window level to better define the boundaries of the structure relative to its surroundings.
- Use Post Processing to smoothen and fill in cavities.
- The following shortcut keys may be helpful in reducing contouring time.
 - '2' to toggle between 2D/3D-brush.
 - 'a' to toggle between static/adaptive brush.
 - Hold the 'shift' key to use brush as an eraser.
 - Press down on the 'mouse wheel' & move the mouse to the left/right to decrease/increase brush diameter.
- Copy and Paste a contour on a slice by switching to a pencil tool.
 - Select Pencil Tool → 'control' key + 'c' → 'control' key + 'v'.
- Delete contour on current plane.
 - Select Pencil Tool → 'delete' key.
- Switch between mouse cursor, pencil, brush, and eraser tool.
 - Highlight Structure → 'right click' on desired tool.

39


Adaptive Brush Contouring Cont.



40

DOSEPEDIA

Adaptive Brush Contouring Cont.




41


Brachial Plexus Axial Slice-by-Slice

DOSEPEDIA


42




43



44




45



46

ASM & MSM Axial Slice-by-Slice DOSEPEDIA

47



48

Thank You
DOSEPEDIA